

Product Datasheet

Recombinant Mouse FGF-2/FGF basic/bFGF Protein (orb1550858)

Description	Recombinant Mouse FGF-2/FGF basic/bFGF Protein
Species/Host	E. coli
Reactivity	Mouse
Endotoxins	< 1 EU/ug by LAL method.
Target	FGF2
Note	For research use only
Purity	> 95% by SDS-PAGE.
Source	KMP1631, Recombinant Mouse FGF-2/FGF basic/bFGF Protein is produced by E. coli expression system. The target protein is expressed with sequence (Met1-Ser154) of mouse FGF-2/FGF basic/bFGF (Accession #P15655).
Expiration Date	6 months from date of receipt.