

Product Datasheet

Recombinant Mouse FGF1/FGFa/FGF acidic Protein (orb1550860)

Description	Recombinant Mouse FGF1/FGFa/FGF acidic Protein
Species/Host	E. coli
Reactivity	Mouse
Endotoxins	< 1 EU/ug by LAL method.
Target	FGF1
Note	For research use only
Purity	> 95% by SDS-PAGE.
Source	KMP1629, Recombinant Mouse FGF1/FGFa/FGF acidic Protein is produced by E. coli expression system. The target protein is expressed with sequence (Phe16-Asp155) of mouse FGF1/FGFa/FGF acidic (Accession #P61148).
Expiration Date	6 months from date of receipt.