

Product Datasheet

Recombinant Mouse IL-1b/IL-1 beta/IL-1F2 Protein (orb1550861)

Description	Recombinant Mouse IL-1b/IL-1 beta/IL-1F2 Protein
Species/Host	E. coli
Reactivity	Mouse
Endotoxins	< 1 EU/ug by LAL method.
Target	IL1b
Note	For research use only
Purity	> 95% by SDS-PAGE.
Source	KMP1628, Recombinant Mouse IL-1b/IL-1 beta/IL-1F2 Protein is produced by E. coli expression system. The target protein is expressed with sequence (Val118-Ser269) of mouse IL-1b/IL-1 beta/IL-1F2 (Accession #P10749).
Expiration Date	6 months from date of receipt.