

## Product Datasheet

# Recombinant Human FABP3/H-FABP Protein, His Tag (orb1550862)

<b>Description</b>	Recombinant Human FABP3/H-FABP Protein, His Tag
<b>Species/Host</b>	E. coli
<b>Reactivity</b>	Human
<b>Endotoxins</b>	< 1 EU/ug by LAL method.
<b>Target</b>	FABP3
<b>Note</b>	For research use only
<b>Purity</b>	> 95% by SDS-PAGE.
<b>Source</b>	KMP1627, Recombinant Human FABP3/H-FABP Protein is produced by E. coli expression system. The target protein is expressed with sequence (Val2-Ala133) of human FABP3/H-FABP (Accession #P05413) fused with an initial Met at the N-terminus and a 6×His Tag at the N-terminus.
<b>Expiration Date</b>	6 months from date of receipt.