

## Product Datasheet

# Recombinant Human RANK/TNFRSF11A/CD265 Protein, His Tag (orb1550871)

|                        |   |
|------------------------|---|
| <b>Description</b>     | Recombinant Human RANK/TNFRSF11A/CD265 Protein, His Tag   |
| <b>Species/Host</b>    | Human   |
| <b>Reactivity</b>      | Human   |
| <b>Endotoxins</b>      | < 1 EU/ug by LAL method.  |
| <b>Target</b>          | RANK  |
| <b>Note</b>            | For research use only   |
| <b>Purity</b>          | > 95% by SDS-PAGE.  |
| <b>Source</b>          | KMP1618, Recombinant Human RANK/TNFRSF11A/CD265 Protein is produced by Human Cell expression system. The target protein is expressed with sequence (Ile30-Pro212) of human RANK/TNFRSF11A/CD265 (Accession #Q9Y6Q6) fused with a 6×His Tag at the C-terminus. |
| <b>Expiration Date</b> | 6 months from date of receipt.  |