

Product Datasheet

Recombinant Human TGF-beta 1/TGFB1 Protein (orb1550873)

Description	Recombinant Human TGF-beta 1/TGFB1 Protein
Species/Host	Human
Reactivity	Human
Endotoxins	< 1 EU/ug by LAL method.
Target	TGFbeta 1
Note	For research use only
Purity	> 95% by SDS-PAGE.
Source	KMP1616, Recombinant Human TGF-beta 1/TGFB1 Protein is produced by Human Cell expression system. The target protein is expressed with sequence (Ala279-Ser390) of human TGF-beta 1/TGFB1 (Accession #P01137).
Expiration Date	6 months from date of receipt.