

## Product Datasheet

# Recombinant Human TGF-beta 2/TGFB2 Protein (orb1550874)

<b>Description</b>	Recombinant Human TGF-beta 2/TGFB2 Protein
<b>Species/Host</b>	Human
<b>Reactivity</b>	Human
<b>Endotoxins</b>	< 1 EU/ug by LAL method.
<b>Target</b>	TGFbeta 2
<b>Note</b>	For research use only
<b>Purity</b>	> 95% by SDS-PAGE.
<b>Source</b>	KMP1615, Recombinant Human TGF-beta 2/TGFB2 Protein is produced by Human Cell expression system. The target protein is expressed with sequence (Ala303-Ser414) of human TGF-beta 2/TGFB2 (Accession #P61812).
<b>Expiration Date</b>	6 months from date of receipt.