

Product Datasheet

Recombinant Human CD112/Nectin-2/Pvrl2/Pvr/Mph/Pvs Protein, His Tag (orb1550876)

Description	Recombinant Human CD112/Nectin-2/Pvrl2/Pvr/Mph/Pvs Protein, His Tag
Species/Host	Human
Reactivity	Human
Endotoxins	< 1 EU/ug by LAL method.
Target	CD112
Note	For research use only
Purity	> 95% by SDS-PAGE.
Source	KMP1613, Recombinant Human CD112/Nectin-2/Pvrl2/Pvr/Mph/Pvs Protein is produced by Human Cell expression system. The target protein is expressed with sequence (Gln32-Leu360) of human CD112/Nectin-2/Pvrl2/Pvr/Mph/Pvs (Accession #Q92692) fused with a 6×His Tag at the C-terminus.
Expiration Date	6 months from date of receipt.