

Product Datasheet

Recombinant Human NKG2D/CD314/KLRK1 Protein, His Tag (orb1550880)

Description	Recombinant Human NKG2D/CD314/KLRK1 Protein, His Tag
Species/Host	Human
Reactivity	Human
Endotoxins	< 1 EU/ug by LAL method.
Target	NKG2D
Note	For research use only
Purity	> 95% by SDS-PAGE.
Source	KMP1609, Recombinant Human NKG2D/CD314/KLRK1 Protein is produced by Human Cell expression system. The target protein is expressed with sequence (Phe78-Val216) of human NKG2D/CD314/KLRK1 (Accession #P26718) fused with a 6×His Tag at the N-terminus.
Expiration Date	6 months from date of receipt.