

Product Datasheet

Recombinant Human MAG/Siglec-4a Protein, His Tag (orb1550881)

| | |
|------------------------|---|
| Description | Recombinant Human MAG/Siglec-4a Protein, His Tag |
| Species/Host | Human |
| Reactivity | Human |
| Endotoxins | < 1 EU/ug by LAL method. |
| Target | MAG |
| Note | For research use only |
| Purity | > 95% by SDS-PAGE. |
| Source | KMP1608, Recombinant Human MAG/Siglec-4a Protein is produced by Human Cell expression system. The target protein is expressed with sequence (Gly20-Pro516) of human MAG/Siglec-4a (Accession #P20916) fused with a 6×His Tag at the C-terminus. |
| Expiration Date | 6 months from date of receipt. |