

Product Datasheet

Recombinant Human Activin A/INHBA Protein (orb1550883)

Description	Recombinant Human Activin A/INHBA Protein
Species/Host	Human
Reactivity	Human
Endotoxins	< 1 EU/ug by LAL method.
Target	Activin A
Note	For research use only
Purity	> 95% by SDS-PAGE.
Source	KMP1606, Recombinant Human Activin A/INHBA Protein is produced by Human Cell expression system. The target protein is expressed with sequence (Gly311-Ser426) of human Activin A/INHBA (Accession #P08476).
Expiration Date	6 months from date of receipt.