

Product Datasheet

Recombinant Human IL-25/IL-17E Protein, His Tag (orb1550884)

Catalog Number	orb1550884
Category	Proteins
Description	Recombinant Human IL-25/IL-17E Protein, His Tag
Target	IL25
Reactivity	Human
Purity	> 95% by SDS-PAGE.
Expression System	HEK293 Cells
Endotoxins	> 95% by SDS-PAGE.
Source	Recombinant Human IL-25/IL-17E Protein is produced by Human Cell expression system. The target protein is expressed with sequence (Tyr33-Gly177) of human IL-25/IL-17E (Accession #Q9H293) fused with a 6His Tag at the C-terminus.
Note	For research use only

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)