

Product Datasheet

Recombinant Human CD155/PVR Protein, His Tag (orb1550906)

Description	Recombinant Human CD155/PVR Protein, His Tag
Species/Host	Human
Reactivity	Human
Endotoxins	< 1.0 EU/ug by LAL method.
Target	CD155
Note	For research use only
Purity	> 97% by SDS-PAGE.
Source	KMP1583, Recombinant Human CD155/PVR Protein is produced by HEK293 Cells expression system. The target protein is expressed with sequence (Asp28-Asn343) of human CD155/PVR (Accession #NP_006496.4) fused with a 6×His Tag at the C-terminus.
Expiration Date	6 months from date of receipt.