

Product Datasheet

Recombinant Human E-Selectin/CD62E Protein, His Tag (orb1550910)

Description	Recombinant Human E-Selectin/CD62E Protein, His Tag
Species/Host	Human
Reactivity	Human
Endotoxins	Please contact us for more information.
Target	ESelectin
Note	For research use only
Purity	> 97% by SDS-PAGE.
Source	KMP1579, Recombinant Human E-Selectin/CD62E Protein is produced by HEK293 expression system. The target protein is expressed with sequence (Trp22-Pro556) of human E-Selectin/CD62E (Accession #NP_000441.2) fused with a 6×His Tag at the C-terminus.
Expiration Date	6 months from date of receipt.