

## Product Datasheet

# Recombinant Human ICAM-2 Protein, Fc & His Tag (orb1550913)

<b>Description</b>	Recombinant Human ICAM-2 Protein, Fc & His Tag
<b>Species/Host</b>	Human
<b>Reactivity</b>	Human
<b>Endotoxins</b>	Please contact us for more information.
<b>Target</b>	ICAM2
<b>Note</b>	For research use only
<b>Purity</b>	> 97% by SDS-PAGE.
<b>Source</b>	KMP1576, Recombinant Human ICAM-2 Protein is produced by HEK293 expression system. The target protein is expressed with sequence (Lys25-Gln216) of human ICAM-2 (Accession #NP_000864.2) fused with an Fc, 6×His Tag at the C-terminus.
<b>Expiration Date</b>	6 months from date of receipt.