

Product Datasheet

Recombinant Human BTN1A1/Butyrophilin Protein, Fc & His Tag (orb1550915)

Description	Recombinant Human BTN1A1/Butyrophilin Protein, Fc & His Tag
Species/Host	Human
Reactivity	Human
Endotoxins	< 1.0 EU/ug by LAL method.
Target	BTN1A1
Note	For research use only
Purity	> 97% by SDS-PAGE.
Source	KMP1574, Recombinant Human BTN1A1/Butyrophilin Protein is produced by HEK293 expression system. The target protein is expressed with sequence (Ala 27 - Arg 242) of human BTN1A1/Butyrophilin (Accession #NP_001723.2) fused with an Fc, 6xHis Tag at the C-terminus.
Expiration Date	6 months from date of receipt.