

Product Datasheet

Recombinant Human CD300C Protein, His Tag (orb1550918)

Description	Recombinant Human CD300C Protein, His Tag
Species/Host	Human
Reactivity	Human
Endotoxins	< 1.0 EU/ug by LAL method.
Target	CD300C
Note	For research use only
Purity	> 95% by SDS-PAGE.
Source	KMP1571, Recombinant Human CD300C Protein is produced by HEK293 Cells expression system. The target protein is expressed with sequence (Gly 21 - Arg 183) of human CD300C (Accession #NP_006669.1) fused with a 6×His Tag at the C-terminus.
Expiration Date	6 months from date of receipt.