

Product Datasheet

Recombinant Human CD46 Protein, His Tag (orb1550920)

Description	Recombinant Human CD46 Protein, His Tag
Species/Host	Human
Reactivity	Human
Endotoxins	< 0.1 EU/ug by LAL method.
Target	CD46
Note	For research use only
Purity	> 90% by SDS-PAGE.
Source	KMP1569, Recombinant Human CD46 Protein is produced by HEK293 Cells expression system. The target protein is expressed with sequence (Cys35- Asp328) of human CD46 (Accession #NP_758861.1) fused with a 6×His Tag at the C-terminus.
Expiration Date	6 months from date of receipt.