

Product Datasheet

Recombinant Human TGFBR2 Protein, Fc & His Tag (orb1550924)

Description	Recombinant Human TGFBR2 Protein, Fc & His Tag
Species/Host	Human
Reactivity	Human
Endotoxins	< 0.1 EU/ug by LAL method.
Target	TGFBR2
Note	For research use only
Purity	> 95% by SDS-PAGE.
Source	KMP1565, Recombinant Human TGFBR2 Protein is produced by HEK293 expression system. The target protein is expressed with sequence (Thr23-Asp159) of human TGFBR2 (Accession #NP_003233.4) fused with an Fc, 6×His Tag at the C-terminus.
Expiration Date	6 months from date of receipt.