

Product Datasheet

Recombinant Human GPNMB Protein, Fc & His Tag (orb1550929)

Description	Recombinant Human GPNMB Protein, Fc & His Tag
Species/Host	Human
Reactivity	Human
Endotoxins	< 0.1 EU/ug by LAL method.
Target	GPNMB
Note	For research use only
Purity	> 90% by SDS-PAGE.
Source	KMP1560, Recombinant Human GPNMB Protein is produced by HEK293 expression system. The target protein is expressed with sequence (Met1-Pro474) of human GPNMB (Accession #NP_002501.1) fused with an Fc, 6×His Tag at the C-terminus.
Expiration Date	6 months from date of receipt.