

Product Datasheet

Recombinant Human TIM-3 Protein, Fc & His Tag (orb1550931)

Description	Recombinant Human TIM-3 Protein, Fc & His Tag
Species/Host	Human
Reactivity	Human
Endotoxins	< 0.1 EU/ug by LAL method.
Target	TIM3
Note	For research use only
Purity	> 87% by SDS-PAGE.
Source	KMP1558, Recombinant Human TIM-3 Protein is produced by HEK293 expression system. The target protein is expressed with sequence (Ser22-Arg200) of human TIM-3 (Accession #NP_116171.3) fused with an Fc, 6×His Tag at the C-terminus.
Expiration Date	6 months from date of receipt.