

Product Datasheet

Recombinant Human CNTF R alpha Protein, His Tag (orb1550935)

Description	Recombinant Human CNTF R alpha Protein, His Tag
Species/Host	Human
Reactivity	Human
Endotoxins	< 0.1 EU/ug by LAL method.
Target	CNTF R alpha
Note	For research use only
Purity	> 95% by SDS-PAGE.
Source	KMP1553, Recombinant Human CNTF R alpha Protein is produced by HEK293 Cells expression system. The target protein is expressed with sequence (Gln 23 - Pro 346) of human CNTF R alpha (Accession #NP_001833.1) fused with a 6×His Tag at the C-terminus.
Expiration Date	6 months from date of receipt.