

Product Datasheet

Recombinant Human BTN1A1/Butyrophilin Protein, His Tag (orb1550942)

Description	Recombinant Human BTN1A1/Butyrophilin Protein, His Tag
Species/Host	Human
Reactivity	Human
Endotoxins	< 0.1 EU/ug by LAL method.
Target	BTN1A1
Note	For research use only
Purity	> 97% by SDS-PAGE.
Source	KMP1546, Recombinant Human BTN1A1/Butyrophilin Protein is produced by HEK293 Cells expression system. The target protein is expressed with sequence (Ala 27 - Arg 242) of human BTN1A1/Butyrophilin (Accession #NP_001723.2) fused with a 6×His Tag at the C-terminus.
Expiration Date	6 months from date of receipt.