

Product Datasheet

Recombinant Human ALK-7 Protein, Fc & His Tag (orb1550944)

Description	Recombinant Human ALK-7 Protein, Fc & His Tag
Species/Host	Human
Reactivity	Human
Endotoxins	< 0.1 EU/ug by LAL method.
Target	ALK7
Note	For research use only
Purity	> 95% by SDS-PAGE.
Source	KMP1544, Recombinant Human ALK-7 Protein is produced by HEK293 expression system. The target protein is expressed with sequence (Leu22-Glu113) of human ALK-7 Protein (Accession #NP_660302.2) fused with an Fc, 6×His Tag at the C-terminus.
Expiration Date	6 months from date of receipt.