

Product Datasheet

Recombinant Human CXADR/CAR Protein, Fc & His Tag (orb1550946)

Description	Recombinant Human CXADR/CAR Protein, Fc & His Tag
Species/Host	Human
Reactivity	Human
Endotoxins	< 0.1 EU/ug by LAL method.
Target	CXADR
Note	For research use only
Purity	> 95% by SDS-PAGE.
Source	KMP1542, Recombinant Human CXADR/CAR Protein is produced by HEK293 expression system. The target protein is expressed with sequence (Leu20 - Gly237) of human CXADR/CAR (Accession #NP_001329.1) fused with an Fc, 6×His Tag at the C-terminus.
Expiration Date	6 months from date of receipt.