

Product Datasheet

Recombinant Human CLEC10A/MGL1/CD301 Protein, His Tag (orb1550973)

Description	Recombinant Human CLEC10A/MGL1/CD301 Protein, His Tag
Species/Host	Human
Reactivity	Human
Endotoxins	< 1.0 EU/ug by LAL method.
Target	CLEC10A
Note	For research use only
Purity	> 90% by SDS-PAGE.
Source	KMP1514, Recombinant Human CLEC10A/MGL1/CD301 Protein is produced by HEK293 Cells expression system. The target protein is expressed with sequence (Gln61-His Tag316) of human CLEC10A/MGL1/CD301 (Accession #NP_878910.1) fused with a 6×His Tag at the C-terminus.
Expiration Date	6 months from date of receipt.