

Product Datasheet

Recombinant Human IGF2 Protein (orb1550979)

Description	Recombinant Human IGF2 Protein
Species/Host	E. coli
Reactivity	Human
Endotoxins	< 1 EU/ug by LAL method.
Target	IGF2
Note	For research use only
Purity	> 95% by SDS-PAGE.
Source	KMP1508, Recombinant Human IGF2 Protein is produced by E. coli expression system. The target protein is expressed with sequence (Ala25-Glu91) of human IGF2 (Accession #P01344).
Expiration Date	6 months from date of receipt.