

Product Datasheet

Recombinant Human FGF-9 Protein (orb1550982)

Description	Recombinant Human FGF-9 Protein
Species/Host	E. coli
Reactivity	Human
Endotoxins	< 1 EU/ug by LAL method.
Target	FGF9
Note	For research use only
Purity	> 95% by SDS-PAGE.
Source	KMP1505, Recombinant Human FGF-9 Protein is produced by E. coli expression system. The target protein is expressed with sequence (Met1-Ser208) of human FGF-9 (Accession #P31371).
Expiration Date	6 months from date of receipt.