

Product Datasheet

Recombinant Human Fibronectin/FN/FN1 Protein (orb1550983)

Description	Recombinant Human Fibronectin/FN/FN1 Protein
Species/Host	E. coli
Reactivity	Human
Endotoxins	< 1 EU/ug by LAL method.
Target	Fibronectin
Note	For research use only
Purity	> 95% by SDS-PAGE.
Source	KMP1504, Recombinant Human Fibronectin/FN/FN1 Protein is produced by E. coli expression system. The target protein is expressed with sequence (Pro1270-Ser1546 & Ala1721-Thr2016) of human Fibronectin/FN/FN1 (Accession #P02751).
Expiration Date	6 months from date of receipt.