

Product Datasheet

Recombinant Human FGF-4 Protein (orb1550985)

Description	Recombinant Human FGF-4 Protein
Species/Host	E. coli
Reactivity	Human
Endotoxins	< 1 EU/ug by LAL method.
Target	FGF4
Note	For research use only
Purity	> 95% by SDS-PAGE.
Source	KMP1502, Recombinant Human FGF-4 Protein is produced by E. coli expression system. The target protein is expressed with sequence (Ser54-Leu206) of human FGF-4 (Accession #P08620) fused with an initial Met at the N-terminus.
Expiration Date	6 months from date of receipt.