

Product Datasheet

Recombinant Human DKK1 Protein, His Tag (orb1550987)

Description	Recombinant Human DKK1 Protein, His Tag
Species/Host	Human
Reactivity	Human
Endotoxins	< 1 EU/ug by LAL method.
Target	DKK1
Note	For research use only
Purity	> 95% by SDS-PAGE.
Source	KMP1500, Recombinant Human DKK1 Protein is produced by Human Cell expression system. The target protein is expressed with sequence (Thr32-His Tag266) of human DKK1 (Accession #O94907) fused with an 8×His Tag at the N-terminus.
Expiration Date	6 months from date of receipt.