

## Product Datasheet

# Recombinant Human Cadherin-1/E-Cadherin/CDH1 Protein, His Tag (orb1550992)

<b>Description</b>	Recombinant Human Cadherin-1/E-Cadherin/CDH1 Protein, His Tag
<b>Species/Host</b>	Human
<b>Reactivity</b>	Human
<b>Endotoxins</b>	< 1 EU/ug by LAL method.
<b>Target</b>	Cadherin1
<b>Note</b>	For research use only
<b>Purity</b>	> 95% by SDS-PAGE.
<b>Source</b>	KMP1495, Recombinant Human Cadherin-1/E-Cadherin/CDH1 Protein is produced by Human Cell expression system. The target protein is expressed with sequence (Asp155-Ile707) of human Cadherin-1/E-Cadherin/CDH1 (Accession #P12830) fused with a 6×His Tag at the C-terminus.
<b>Expiration Date</b>	6 months from date of receipt.