

Product Datasheet

Recombinant Human CCL5/RANTES Protein (orb1550993)

Description	Recombinant Human CCL5/RANTES Protein
Species/Host	E. coli
Reactivity	Human
Endotoxins	< 1 EU/ug by LAL method.
Target	CCL5
Note	For research use only
Purity	> 95% by SDS-PAGE.
Source	KMP1494, Recombinant Human CCL5/RANTES Protein is produced by E. coli expression system. The target protein is expressed with sequence (Ser24-Ser91) of human CCL5/RANTES (Accession #P13501).
Expiration Date	6 months from date of receipt.