

Product Datasheet

Recombinant Human CCL3/MIP-1 alpha Protein (orb1550994)

Description	Recombinant Human CCL3/MIP-1 alpha Protein
Species/Host	E. coli
Reactivity	Human
Endotoxins	< 1 EU/ug by LAL method.
Target	CCL3
Note	For research use only
Purity	> 95% by SDS-PAGE.
Source	KMP1493, Recombinant Human CCL3/MIP-1 alpha Protein is produced by E. coli expression system. The target protein is expressed with sequence (Ser24-Ala92) of human CCL3/MIP-1 alpha (Accession #P10147).
Expiration Date	6 months from date of receipt.