

Product Datasheet

Recombinant Human CD27/TNFRSF7 Protein, Fc & His Tag (orb1550995)

Description	Recombinant Human CD27/TNFRSF7 Protein, Fc & His Tag
Species/Host	Human
Reactivity	Human
Endotoxins	< 1 EU/ug by LAL method.
Target	CD27
Note	For research use only
Purity	> 95% by SDS-PAGE.
Source	KMP1492, Recombinant Human CD27/TNFRSF7 Protein is produced by Human Cell expression system. The target protein is expressed with sequence (Thr21-Ile192) of human CD27/TNFRSF7 (Accession #P26842) fused with an Fc, 6×His Tag at the C-terminus.
Expiration Date	6 months from date of receipt.