

Product Datasheet

Recombinant Human BDNF Protein (orb1550999)

Description	Recombinant Human BDNF Protein
Species/Host	E. coli
Reactivity	Human
Endotoxins	< 1 EU/ug by LAL method.
Target	BDNF
Note	For research use only
Purity	> 95% by SDS-PAGE.
Source	KMP1488, Recombinant Human BDNF Protein is produced by E. coli expression system. The target protein is expressed with sequence (His Tag129-Arg247) of human BDNF (Accession #P23560) fused with an initial Met at the N-terminus.
Expiration Date	6 months from date of receipt.