

## Product Datasheet

### Recombinant Human LGALS7 Protein (orb1551005)

<b>Description</b>	Recombinant Human LGALS7 Protein
<b>Species/Host</b>	E. coli
<b>Reactivity</b>	Human
<b>Endotoxins</b>	< 1 EU/ug by LAL method.
<b>Target</b>	LGALS7
<b>Note</b>	For research use only
<b>Purity</b>	> 95% by SDS-PAGE.
<b>Source</b>	KMP1482, Recombinant Human LGALS7 Protein is produced by E.coli expression system. The target protein is expressed with sequence (Met1-Phe136) of human LGALS7 (Accession #P47929).
<b>Expiration Date</b>	6 months from date of receipt.