

## Product Datasheet

# Recombinant Human IGSF11 Protein, Fc His Tag (orb1551006)

<b>Description</b>	Recombinant Human IGSF11 Protein, Fc His Tag
<b>Species/Host</b>	Mammal
<b>Reactivity</b>	Human
<b>Endotoxins</b>	< 1 EU/ug by LAL method.
<b>Target</b>	IGSF11
<b>Note</b>	For research use only
<b>Purity</b>	> 95% by SDS-PAGE.
<b>Source</b>	KMP1481, Recombinant Human IGSF11 Protein is produced by Mammalian expression system. The target protein is expressed with sequence (Leu23-Gly245) of human IGSF11 (Accession #Q5DX21) fused with an Fc tag at the C-terminus.
<b>Expiration Date</b>	6 months from date of receipt.