

Product Datasheet

Recombinant Human CD72 Protein, Trx & His Tag (orb1551007)

Description	Recombinant Human CD72 Protein, Trx & His Tag
Species/Host	E. coli
Reactivity	Human
Endotoxins	< 1 EU/ug by LAL method.
Target	CD72
Note	For research use only
Purity	> 90% by SDS-PAGE.
Source	KMP1480, Recombinant Human CD72 Protein is produced by E.coli expression system. The target protein is expressed with sequence (Arg117-Cys226) of human CD72 (Accession #P21854) fused with a Trx, 6xHis Tag at the N-terminus.
Expiration Date	6 months from date of receipt.