

Product Datasheet

Recombinant Human IGF2BP-2 Protein, T7 & His Tag (orb1551008)

Description	Recombinant Human IGF2BP-2 Protein, T7 & His Tag
Species/Host	E. coli
Reactivity	Human
Endotoxins	< 1 EU/ug by LAL method.
Target	IGF2BP2
Note	For research use only
Purity	> 80% by SDS-PAGE.
Source	KMP1479, Recombinant Human IGF2BP-2 Protein is produced by E.coli expression system. The target protein is expressed with sequence (Met1-Thr220) of human IGF2BP-2 (Accession #Q9Y6M1) fused with a T7 tag at the N-terminus and a 6xHis Tag at the C-terminus.
Expiration Date	6 months from date of receipt.