

Product Datasheet

Recombinant Human S100A13 Protein (orb1551015)

Description	Recombinant Human S100A13 Protein
Species/Host	E. coli
Reactivity	Human
Endotoxins	< 1 EU/ug by LAL method.
Target	S100A13
Note	For research use only
Purity	> 95% by SDS-PAGE.
Source	KMP1472, Recombinant Human S100A13 Protein is produced by E.coli expression system. The target protein is expressed with sequence (Ala2-Lys98) of human S100A13 (Accession #Q99584).
Expiration Date	6 months from date of receipt.