

## Product Datasheet

# Recombinant Human TXN2 Protein, His Tag (orb1551016)

<b>Description</b>	Recombinant Human TXN2 Protein, His Tag
<b>Species/Host</b>	E. coli
<b>Reactivity</b>	Human
<b>Endotoxins</b>	< 1 EU/ug by LAL method.
<b>Target</b>	TXN2
<b>Note</b>	For research use only
<b>Purity</b>	> 95% by SDS-PAGE.
<b>Source</b>	KMP1471, Recombinant Human TXN2 Protein is produced by E.coli expression system. The target protein is expressed with sequence (Thr60-Gly166) of human TXN2 (Accession #Q99757) fused with a 6xHis Tag at the N-terminus.
<b>Expiration Date</b>	6 months from date of receipt.