

Product Datasheet

Recombinant Human SUMO1 Protein, His Tag (orb1551017)

Description	Recombinant Human SUMO1 Protein, His Tag
Species/Host	E. coli
Reactivity	Human
Endotoxins	< 1 EU/ug by LAL method.
Target	SUMO1
Note	For research use only
Purity	> 95% by SDS-PAGE.
Source	KMP1470, Recombinant Human SUMO1 Protein is produced by E.coli expression system. The target protein is expressed with sequence (Met1-Val101) of human SUMO1 (Accession #AAH66306) fused with a 6xHis Tag at the N-terminus.
Expiration Date	6 months from date of receipt.