

Product Datasheet

Recombinant Human FABP2 Protein, His Tag (orb1551019)

Description	Recombinant Human FABP2 Protein, His Tag
Species/Host	E. coli
Reactivity	Human
Endotoxins	< 1 EU/ug by LAL method.
Target	FABP2
Note	For research use only
Purity	> 95% by SDS-PAGE.
Source	KMP1468, Recombinant Human FABP2 Protein is produced by E.coli expression system. The target protein is expressed with sequence (Met1-Asp132) of human FABP2 (Accession #P12104) fused with a 6xHis Tag at the N-terminus and a 6xHis Tag at the C- terminus.
Expiration Date	6 months from date of receipt.