

## Product Datasheet

# Recombinant Human LAIR1 Protein, His Tag (orb1551023)

<b>Description</b>	Recombinant Human LAIR1 Protein, His Tag
<b>Species/Host</b>	Mammal
<b>Reactivity</b>	Human
<b>Endotoxins</b>	< 1 EU/ug by LAL method.
<b>Target</b>	LAIR1
<b>Note</b>	For research use only
<b>Purity</b>	> 95% by SDS-PAGE.
<b>Source</b>	KMP1464, Recombinant Human LAIR1 Protein is produced by Mammalian expression system. The target protein is expressed with sequence (Gln22-His Tag163) of human LAIR1 (Accession #Q6GTX8) fused with a 6xHis Tag at the C-terminus.
<b>Expiration Date</b>	6 months from date of receipt.