

Product Datasheet

Recombinant Human Noggin / NOG Protein, Fc His Tag (orb1551025)

Description	Recombinant Human Noggin / NOG Protein, Fc His Tag
Species/Host	Human
Reactivity	Human
Endotoxins	< 1.0 EU/ug by LAL method.
Target	Noggin
Note	For research use only
Purity	> 95 % by SDS-PAGE.
Source	KMP1462, Recombinant Human Noggin / NOG Protein is produced by HEK293 Cells expression system. The target protein is expressed with sequence (Met1-Cys232) of human Noggin / NOG (Accession #NP_005441.1) fused with an Fc tag at the C-terminus.
Expiration Date	6 months from date of receipt.