

## Product Datasheet

### Recombinant Human CXCL4 Protein (orb1551026)

|                        |  |
|------------------------|--|
| <b>Description</b>     | Recombinant Human CXCL4 Protein  |
| <b>Species/Host</b>    | E. coli  |
| <b>Reactivity</b>      | Human  |
| <b>Endotoxins</b>      | Please contact us for more information.  |
| <b>Target</b>          | CXCL4  |
| <b>Note</b>            | For research use only  |
| <b>Purity</b>          | > 80 % by SDS-PAGE.  |
| <b>Source</b>          | KMP1461, Recombinant Human CXCL4 Protein is produced by E. coli expression system. The target protein is expressed with sequence (Glu32-Ser101) of human CXCL4 (Accession #NP_002610.1). |
| <b>Expiration Date</b> | 6 months from date of receipt.   |