

Product Datasheet

Recombinant Human VSIG8 Protein, His Tag (orb1551027)

Description	Recombinant Human VSIG8 Protein, His Tag
Species/Host	Mammal
Reactivity	Human
Endotoxins	< 1 EU/ug by LAL method.
Target	VSIG8
Note	For research use only
Purity	> 95% by SDS-PAGE.
Source	KMP1460, Recombinant Human VSIG8 Protein is produced by Mammalian expression system. The target protein is expressed with sequence (Val22-Gly263) of human VSIG8 (Accession #Q5VU13) fused with a 6xHis Tag at the C- terminus.
Expiration Date	6 months from date of receipt.