

Product Datasheet

Recombinant Mouse TNFR2 / TNFRSF1B Protein, Fc His Tag (orb1551028)

Description	Recombinant Mouse TNFR2 / TNFRSF1B Protein, Fc His Tag
Species/Host	Human
Reactivity	Mouse
Endotoxins	< 1.0 EU/ug by LAL method.
Target	TNFR2
Note	For research use only
Purity	> 95 % by SDS-PAGE.
Source	KMP1459, Recombinant Mouse TNFR2 / TNFRSF1B Protein is produced by HEK293 Cells expression system. The target protein is expressed with sequence (Met1-Gly258) of mouse TNFR2 / TNFRSF1B (Accession #NP_035740.2) fused with an Fc tag at the C-terminus.
Expiration Date	6 months from date of receipt.