

Product Datasheet

Recombinant Human Insulin Receptor Protein, His Tag (orb1551030)

Description	Recombinant Human Insulin Receptor Protein, His Tag
Species/Host	Human
Reactivity	Human
Endotoxins	< 1.0 EU/ug by LAL method.
Target	Insulin Receptor
Note	For research use only
Purity	> 95 % by SDS-PAGE.
Source	KMP1457, Recombinant Human Insulin Receptor Protein is produced by HEK293 Cells expression system. The target protein is expressed with sequence (Met1-Lys956) of human Insulin Receptor (Accession #NP_000199.2) fused with a 6×His Tag at the C-terminus.
Expiration Date	6 months from date of receipt.